

## Weekday Manufactured

Expected mCi of Tc 99m at 8:00 a.m. ET

| Gen Size | 1 day Post Cal | 2 day Post Cal | 3 day Post Cal | 4 day Post Cal | 5 day Post Cal | 6 day Post Cal | 7 day Post Cal | 8 day Post Cal | 9 day Post Cal | 10 day Post Cal | 11 day Post Cal | 12 day Post Cal | 13 day Post Cal | 14 day Post Cal |
|----------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| 1000     | 610            | 535            | 420            | 326            | 254            | 197            | 153            | 119            | 92             | 72              | 56              | 43              | 34              | 26              |
| 2000     | 1220           | 1070           | 839            | 653            | 507            | 394            | 306            | 238            | 185            | 144             | 112             | 87              | 67              | 52              |
| 2500     | 1525           | 1338           | 1049           | 816            | 634            | 493            | 383            | 297            | 231            | 180             | 139             | 108             | 84              | 65              |
| 3000     | 1830           | 1605           | 1259           | 979            | 761            | 591            | 459            | 357            | 277            | 215             | 167             | 130             | 101             | 79              |
| 4000     | 2440           | 2140           | 1678           | 1305           | 1014           | 788            | 612            | 476            | 370            | 287             | 223             | 173             | 135             | 105             |
| 4500     | 2745           | 2408           | 1888           | 1468           | 1141           | 887            | 689            | 535            | 416            | 323             | 251             | 195             | 152             | 118             |
| 5000     | 3050           | 2675           | 2098           | 1631           | 1268           | 985            | 765            | 595            | 462            | 359             | 279             | 217             | 168             | 131             |
| 6000     | 3660           | 3211           | 2518           | 1958           | 1521           | 1182           | 918            | 714            | 554            | 431             | 335             | 260             | 202             | 157             |
| 7500     | 4575           | 4013           | 3147           | 2447           | 1902           | 1478           | 1148           | 892            | 693            | 539             | 418             | 325             | 253             | 196             |
| 10000    | 6100           | 5351           | 4196           | 3263           | 2535           | 1970           | 1531           | 1189           | 924            | 718             | 558             | 434             | 337             | 262             |
| 12500    | 7626           | 6689           | 5245           | 4078           | 3169           | 2463           | 1913           | 1487           | 1155           | 898             | 697             | 542             | 421             | 327             |
| 15000    | 9151           | 8026           | 6294           | 4894           | 3803           | 2955           | 2296           | 1784           | 1386           | 1077            | 837             | 650             | 505             | 393             |
| 18000    | 10981          | 9632           | 7553           | 5873           | 4564           | 3546           | 2755           | 2141           | 1663           | 1293            | 1004            | 78              | 606             | 471             |

- Calibration is 12 Noon ET on day of manufacture
- Manufacture Day: Monday, Tuesday, and Friday (subject to change)
- Expected mCi of Tc 99m assumes 90% elution efficiency on that day

Please read the accompanying full [Prescribing Information](#) also available at [LantheusSPECT.com](#).